

## UltraMicroscope Specialist Training

# How to clean a cuvette

038\_002\_CuvetteCleaning\_LH Version 1.0

### Materials needed:

- Dish soap
- Water
- Ethanol 70%–99%
- Lint free tissue
- Cuvette (procedure is the same for UltraMicroscope II and UltraMicroscope Blaze™ cuvettes)

1. Soak the cuvette in warm soap water. Do not use any tools like sponges or cloths.



2. Rinse the cuvette with water.



3. Add a small amount of ethanol to remove any oil residues that are left. Carefully turn the cuvette to disperse the ethanol everywhere.







4. Use a lint free tissue with some ethanol to clean the windows on the sides. Make sure to dry of the ethanol, otherwise the windows can get streaky.



UltraMicroscope Blaze™ is a trademark of Miltenyi Biotec GmbH. Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. Copyright © Miltenyi Biotec 2021. All rights reserved.



5. Dry the rest of the cuvette with a paper towel.